

Biological Indicator

H₂O₂ Sterilization Consumables

Biological indicators are used for biological monitoring of Hydrogen Peroxide Low Temperature Plasma Sterilizer. The whole process goes as follows: sterilize the biological indicator together with sterilization articles (please do not break biological indicator when sterilizing) and take out the indicator after sterilization and break it; and then put it into incubator for monitoring spores growth. If with growth, the culture medium color would turn yellow and if not, no color changes on culture medium. The final incubation result will out in corresponding hours.



Yellow color means   **POSITIVE**

FAILURE (Positive: Purple → Yellow)

Purple color means   **NEGATIVE**

SUCCESS (Negative: Purple → Purple)



Model	Description	Package
FUGBI0102-2	Rapid readout BI Incubation time: 20 Minutes	50 Pieces/Box