

# Biological Indicator

## Steam Sterilization Consumables

Ugaiya Steam Biological indicators are a special type of living microbial products that can be used to confirm the performance of sterilization equipment, the verification of sterilization equipment, the verification of sterilization procedures, and the monitoring of sterilization effects in production processes. It is defined as a defined resistance to a specific sterilization process and is available in the inner packaging for the carrier.



**Yellow color means**   **POSITIVE**

**FAILURE** (Positive: Purple → Yellow)

**Purple color means**   **NEGATIVE**

**SUCCESS** (Negative: Purple → Purple)

| Model       | Description                 | Package       |
|-------------|-----------------------------|---------------|
| FUGBI0201-2 | Incubation time: 20 Minutes | 50 Pieces/Box |

